


**Thermocouple magnétique bouton de traction JIS  
Type K 1,9 kg, cble isol PFA avec sur-tresse en  
acier inoxydable se terminant dans une prise  
standard - 2 m**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Spcifications****Specifications**

<b>Product Code</b>	XE-0627-001
<b>General Description</b>	1.9kg Pull Button Magnet Thermocouple
<b>Thermocouple Type</b>	K JIS
<b>Magnet</b>	19.1mm Diameter / 12.7mm height - Max pull: 1.9 KG
<b>Cable Length</b>	2m
<b>Cable Type</b>	PFA Insulated Cable with stainless steel overbraid 7/0.2mm
<b>Termination Type</b>	Standard Plug
<b>Max. Temperature</b>	250°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	JIS (Japanese Industrial Standards)