

**Thermocouple pince crocodile avec cble surtress
en acier inoxydable en fibre de verre - JIS Type K
- Longueur 3m**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-0631-001
General Description	Crocodile Clip Thermocouple with Glassfibre Stainless Steel Overbraided Cable
Thermocouple Type	K JIS
Cable Length	3m
Cable Type	Stainless steel over-braided PFA Teflon®
Cable Termination	Miniature Plug
Max. Temperature	250°C
Min. Temperature	-50°C
Standards Met	JIS (Japanese Industrial Standards)