



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

**Thermocouple mis la terre fabriqué dans un tube en acier inoxydable, câble surstress en acier inoxydable en fibre de verre - Type K JIS - 4,0 mm de diamètre x 13 mm de longueur**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Specifications****Specifications**

<b>Product Code</b>	XE-0646-001
<b>Sensor Type</b>	316 stainless steel tube in various sizes
<b>Thermocouple Type</b>	K JIS
<b>Thermocouple Junction</b>	Grounded at tip for fast response
<b>Cable Length</b>	2m
<b>Cable Type</b>	Glassfibre stainless steel overbraided cable
<b>Cable Termination</b>	Bare Wire Tails
<b>Max. Temperature</b>	350°C
<b>Min. Temperature</b>	-60°C
<b>Standards Met</b>	JIS (Japanese Industrial Standards)
<b>Dimensions</b>	4.0mm Diameter x 13mm Length