




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Thermocouple de surface de tuyau, surtresse en acier inoxydable en fibre de verre - JIS Type K - pour sadapter un diamtre de tuyau de 22 30 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications

Specifications

Product Code	XE-0659-001
General Description	Stainless steel pipe clamp sensor with adjustable diameter
Thermocouple Type	K JIS
Thermocouple Junction	Grounded tip for fast response
Cable Length	2m
Cable Type	Glassfibre stainless steel braided extension cable
Core / Strands	7/0.2mm
Termination Type	Bare Tails
Max. Temperature	350°C
Standards Met	JIS (Japanese Industrial Standards)
To Fit Pipe Diameters	22-30mm