

**ANSI Type K 4.5kg Pull Power (Horseshoe)  
Magnet Thermocouple, PFA Insulated Cable avec  
tresse en acier inoxydable se terminant en  
miniature ou standard**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Specifications****Specifications**

<b>General Description</b>	4.5kg Pull Power (Horseshoe) Magnet Thermocouple
<b>Sensor Type</b>	Magnet Thermocouple
<b>Thermocouple Type</b>	K
<b>Magnet</b>	30 x 20 x 20mm (L x H x W). Max Pull: 4.5kg
<b>Cable Length</b>	1m or 2m
<b>Cable Type</b>	PFA Teflon® insulated with stainless steel over-braid
<b>Core / Strands</b>	7/0.2mm
<b>Termination Type</b>	Miniature or Standard Plug
<b>Max. Temperature</b>	+250°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	ANSI MC96.1