



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

ANSI Type K 9kg Pull Power (Horseshoe) Thermocouple magnétique, cble isol PFA avec tresse en acier inoxydable se terminant en miniature ou fiche standard



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

General Description	Rapid, reliable temperature sensing on ferrous surfaces
Thermocouple Type	K
Thermocouple Junction	Grounded at tip
Magnet	25 x 25 x 40mm (HxWxL) with powerful 9kg pull
Cable Termination	1.5 or 3 metres PFA Teflon® insulated with stainless steel over-braid, fitted plug termination
Termination Type	Fitted Mini or Standard Plug
Max. Temperature	+250°C
Min. Temperature	-50°C
Standards Met	ANSI MC96.1