



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

JIS Type K 4.5kg Pull Power (fer cheval) Thermocouple magnétique, cble isol PFA avec tresse en acier inoxydable se terminant en miniature ou standard



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Spcifications**Specifications**

General Description	4.5kg Pull Power (Horseshoe) Magnet Thermocouple
Sensor Type	Magnet Thermocouple
Thermocouple Type	K
Magnet	30 x 20 x 20mm (L x H x W). Max Pull: 4.5kg
Cable Length	1.0m or 2.0m
Cable Type	Glassfibre insulated lead with stainless steel overbraid
Core / Strands	7/0.2mm
Termination Type	Miniature or Standard Plug
Max. Temperature	+250°C
Min. Temperature	-50°C
Standards Met	JIS (Japanese Industrial Standards)