

**Thermocouple magnétique JIS Type K 9kg Pull Power (fer cheval), cble isol PFA avec sur-tresse en acier inoxydable se terminant par une prise**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Specifications****Specifications**

<b>Product Code</b>	XE-0679-001
<b>General Description</b>	Rapid, reliable temperature sensing on ferrous surfaces
<b>Sensor Type</b>	K JIS
<b>Thermocouple Junction</b>	Grounded at Tip
<b>Magnet</b>	25 x 25 x 40mm (HxWxL) with powerful 9kg pull
<b>Cable Length</b>	3m
<b>Cable Type</b>	PFA Teflon® insulated with stainless steel over-braid, fitted plug termination
<b>Termination Type</b>	Miniature Plug
<b>Max. Temperature</b>	250°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	JIS (Japanese Industrial Standards)