


Standard Rectangulaire Thermocouple Connecteur Fascia Socket IS-K-FF Type K IEC



Standard Fascia Socket (with nylon mounting clip)
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Standard Rectangulaire Thermocouple Connecteur Fascia Socket IS-K-FF Type K IEC

Clip en nylon pour une fixation simple et rapide sur le panneau. Toutes les prises de montage de carnage sont polarises au moyen de broches de taille ingale

Spcifications**Specifications**

Product Code	XE-1407-001
General Description	Accepts standard round pin plug connector. Nylon clip for quick, simple fixing to panel
Thermocouple Type	K IEC
Connector Type	Fascia Mounting Socket
Connector Size	Standard
Colour	Green
Max. Temperature	120°C
+Positive Contact	Nickel Chromium
-Negative Contact	Nickel Alloy
Standards Met	BS EN 50212:1997. RoHS