


# Standard Rectangulaire Thermocouple Connecteur Fascia Socket IS-K-FF Type K IEC



**Standard Fascia Socket (with nylon mounting clip)**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Standard Rectangulaire Thermocouple Connecteur Fascia Socket IS-K-FF Type K IEC**

Clip en nylon pour une fixation simple et rapide sur le panneau. Toutes les prises de montage de carnage sont polarises au moyen de broches de taille ingale

**Spcifications****Specifications**

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | XE-1407-001   |
| <b>General Description</b> | Accepts standard round pin plug connector. Nylon clip for quick, simple fixing to panel |
| <b>Thermocouple Type</b>   | K IEC   |
| <b>Connector Type</b>      | Fascia Mounting Socket  |
| <b>Connector Size</b>      | Standard  |
| <b>Colour</b>              | Green   |
| <b>Max. Temperature</b>    | 120°C   |
| <b>+Positive Contact</b>   | Nickel Chromium   |
| <b>-Negative Contact</b>   | Nickel Alloy  |
| <b>Standards Met</b>       | BS EN 50212:1997. RoHS  |