



# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

**Thermocouple lames ANSI, cble isol en fibre de verre avec sur-tresse en acier inoxydable se terminant par des queues nues, prise miniature**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Specifications

**Specifications**

<b>General Description</b>	Thermocouple mounted on a stainless steel shim
<b>Thermocouple Type</b>	K
<b>Thermocouple Junction</b>	Junction grounded at tip for fast response
<b>Cable Length</b>	1, 2, 3 or 5m
<b>Cable Type</b>	Glassfibre insulated lead with stainless steel overbraid
<b>Core / Strands</b>	7/0.2mm
<b>Termination Type</b>	Bare tails, Miniature or Standard Plug (with cable clamps)
<b>Max. Temperature</b>	+350°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	ANSI MC96.1
<b>Dimensions</b>	25x13mm shim (0.2mm thick)