


**Thermocouple lames, cble isol en fibre de verre
avec sur-tresse en acier inoxydable - Type K JIS
avec terminaison de fiche standard - longueur 1**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Spcifications**Specifications**

Product Code	XE-0698-001
General Description	Thermocouple mounted on a stainless steel shim
Thermocouple Type	K
Thermocouple Junction	Junction grounded at tip for fast response
Cable Length	1M
Cable Type	Glassfibre insulated lead with stainless steel overbraid
Core / Strands	7/0.2mm
Termination Type	Standard Plug (with cable clamps)
Max. Temperature	+350°C
Min. Temperature	-60°C
Standards Met	JIS (Japanese Industrial Standards)
Dimensions	25x13mm shim (0.2mm thick)