



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Thermocouple lames JIS, cble isol en fibre de verre avec sur-tresse en acier inoxydable se terminant par des queues nues, prise miniature



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

General Description	Thermocouple mounted on a stainless steel shim
Thermocouple Type	K
Thermocouple Junction	Junction grounded at tip for fast response
Cable Length	1,2,3 or 5m
Cable Type	Glassfibre insulated lead with stainless steel overbraid
Core / Strands	7/0.2mm
Termination Type	Bare tails, Miniature or Standard Plug (with cable clamps)
Max. Temperature	+350°C
Min. Temperature	-60°C
Standards Met	JIS (Japanese Industrial Standards)
Dimensions	25x13mm shim (0.2mm thick)