




# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

**Sonde thermocouple usage gnral JIS, cble isol  
en fibre de verre avec sur-tresse en acier  
inoxydable se terminant par des queues nues,  
fiche miniature ou standard Type K**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Specifications****Specifications**

<b>General Description</b>	Rugged thermocouple with stainless steel sheath
<b>Thermocouple Type</b>	K
<b>Thermocouple Junction</b>	Grounded hot junction for fast response
<b>Probe Diameter</b>	4.5mm
<b>Probe Length</b>	150mm
<b>Cable Length</b>	1 or 2m
<b>Cable Type</b>	Glassfibre insulated lead with stainless steel overbraid
<b>Core / Strands</b>	7/0.2mm
<b>Termination Type</b>	Bare tails, Miniature or Standard Plug (with cable clamps)
<b>Max. Temperature</b>	+350°C
<b>Min. Temperature</b>	-60°C
<b>Standards Met</b>	JIS