



## Connecteur thermocouple rectangulaire standard Prise de fascia FSTC-CU-FF Type Cu



**Standard Fascia Socket (with nylon mounting clip)**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Standard Rectangular Thermocouple Connector Fascia Socket FSTC-CU-FF Type CU**  
Nylon clip pour une fixation rapide et simple au panneau. Toutes les douilles de montage de fascia sont polarisées au moyen de broches de taille ingale

Spécifications

## Specifications

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | XE-1677-001   |
| <b>General Description</b> | Accepts standard round pin plug connector. Nylon clip for quick, simple fixing to panel |
| <b>Thermocouple Type</b>   | IEC, ANSI, JIS Type CU  |
| <b>Connector Type</b>      | Fascia Mounting Socket  |
| <b>Connector Size</b>      | Standard  |
| <b>Colour</b>              | White   |
| <b>Max. Temperature</b>    | 120°C   |
| <b>+Positive Contact</b>   | Copper  |
| <b>-Negative Contact</b>   | Copper  |
| <b>Standards Met</b>       | BS EN 50212:1997. RoHS  |