

PTFE Olives



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Olives PTFE. Disponible pour sadapter aux diamtres de sonde allant de 1,5 mm 8,0 mm

Specifications

| | |
|----------------------------|--|
| General Description | Olives |
| Material | PTFE |
| Probe Diameter | Suitable for probe diameters between 1.5mm & 8.0mm |