

Sonde pointe plate - Type RTD / PRT - Longueur 2M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Sonde Pt100 usage gnral

Aplati la pointe pour une meilleure rponse thermique

Configuration 2 fils

Classe B, selon IEC 751

Gaine en acier inoxydable 316, 4,0 mm de diamtre x 150 mm de long

Plomb torsad 2 noyaux isol en tflon de 1 mtre de long et isol au tflon

Comprend un raccord compression mle en laiton 1/8" BSPP avec contre-crou

Plage de temprature -50C + 250C

Spcifications

Specifications

Product Code	XE-3698-001
General Description	316 stainless steel probe, 4.0mm diameter (6mm max. across flat area of tip) x 150mm long overall
Sensor Type	Pt100 (100 Ohms @ 0°C) to IEC 751, Class B
Probe Temperature Range	-50°C to +250°C
Cable Length	2M
Cable Type	7/0.2mm Teflon® insulated 2 core twisted leads
Termination Type	2-wire tails
Accessories	1/8" BSPP brass male compression fitting with locknut included