

Sonde bout plat Type RTD / PRT



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Sonde Pt100 usage gnral

Aplati lextmit pour une meilleure rponse thermique

Configuration 2 fils

Classe B, selon CEI 751

Gaine en acier inoxydable 316, diamtre 4,0 mm x 150 mm de long

Plomb torsad isol en tflon de 1 mtre de long 7/0,2 mm 2 conducteurs

Comprend un raccord de compression mle en laiton BSPP de 1/8 po avec crou de verrouillage

Plage de temprature -50C + 250C

Spfications

Specifications

General Description	316 stainless steel probe, 4.0mm diameter (6mm max. across flat area of tip) x 150mm long overall
Sensor Type	Pt100 (100 Ohms @ 0°C) to IEC 751, Class B
Probe Temperature Range	-50°C to +250°C
Cable Length	1m, 2m, 5m
Cable Type	7/0.2mm Teflon® insulated 2 core twisted leads
Termination Type	2-wire tails
Accessories	1/8" BSPP brass male compression fitting with locknut included