


## Terminaux bche thermocouple Type K



Thermocouple Spade Terminal

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

<ul><li>Terminaux bche sertissage</li><li>Prix par paquet de 20 paires </li><li>Convient aux types dtalonnage IEC, ANSI et JIS.</li><li>Fabriqu partir de matriau thermocouple compatible</li><li>Type de thermocouple estampill sur terminal</li><li>Suits Labfacility & autres bornes barrires </li><li>Pour les tailles de fil de 0,5 1,0 mm</li></ul><p><img title="Bande terminale alt="Bande terminale src="{{media url=wysiwyg/barrier\_terminal\_strip.png}} /></p>

## Specifications

<b>Product Code</b>	XE-1982-020
<b>General Description</b>	pack of 20 pairs
<b>+Positive Contact</b>	Nickel Chromium
<b>-Negative Contact</b>	Nickel Alloy