



## Terminales de bÃ¢che thermocouple Type R/S



Thermocouple Spade Terminal

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Bornes a beche sertie

Prix par paquet de 20 paires

Convient aux types d'étalonnage IEC, ANSI et JIS.

Fabriqué a partir d'un matériau thermocouple compatible

Type de thermocouple estampillé sur le terminal

Convient Labfacility et autres terminaux de barriere

Pour les tailles de fil de 0,5 a 1,0 mm

## Specifications

<b>Product Code</b>	XE-1984-020
<b>General Description</b>	pack of 20 pairs
<b>+Positive Contact</b>	Copper
<b>-Negative Contact</b>	Copper Alloy