

## Sonde thermistance usage gnral avec fil PTFE



[Click here to view product](#)

<b>Product Code</b>	XE-9840
<b>General Description</b>	2-wire 10K NTC thermistor housed in a stainless steel sheath
<b>Thermistor</b>	10K NTC Thermistor
<b>BETA 25/85</b>	3977K
<b>Probe Diameter</b>	3mm, 4mm, 4.5mm & 6mm
<b>Cable Length</b>	2M
<b>Cable Type</b>	PTFE insulated twisted leads
<b>Core / Strands</b>	7/0.2mm
<b>Cable Termination</b>	Bare Tails
<b>Max. Temperature</b>	+150Å°C

<b>Min. Temperature</b>	-50Å°C
<b>Accuracy</b>	+/- 0.2K (over 0-70Å°C)
<b>Length</b>	25mm, 50mm, 100mm & 150mm
<b>Number of Wires</b>	2