


Thermistance Capteurs de temprature extrieurs / entrepts frigorifiques



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Thermistance NTC 10K 2 fils dans un boîtier de montage mural résistant aux intempéries IP67 pour une détection précise de la température de l'air extérieur

Caractéristiques principales :

Valeur de résistance : 10K ohms

Type : Thermistance NTC (Negative Temperature Coefficient), ce qui signifie que sa résistance diminue mesure que la température augmente.

Plage de température : Capable de mesurer avec précision des températures de -50C +150C.

Câblage : configuration 2 fils pour une connectivité simple.

Indice IP : IP67

Applications:**Systemes CVC**

: Surveille la température de l'air pour réguler efficacement les systèmes de chauffage, de ventilation et de climatisation.

Spécifications**Specifications**

Product Code	XE-9867-001
General Description	Weatherproof IP67 rated wall mounting Thermistor for accurate outdoor air temperature sensing
Sensor Type	Single Element
Thermistor	10K NTC Thermistor
BETA 25/85	3977K
Sheath Material	Perforated 316 stainless steel probe, 12.7mm x 75mm
Core / Strands	7/0.2mm
Cable Termination	Bare Tails



Gland	M20 cable entry gland
Max. Temperature	+150°C
Min. Temperature	-50°C
IP Rating	67
Accuracy	+/- 0.2K (over 0-70°C)
Number of Wires	2