


Thermistance Capteurs de temprature extrieurs / entrepts frigorifiques



[Click here to view product](#)

Product Code	XE-9867-001
General Description	Weatherproof IP67 rated wall mounting Thermistor for accurate outdoor air temperature sensing
Sensor Type	Single Element
Thermistor	10K NTC Thermistor
BETA 25/85	3977K
Sheath Material	Perforated 316 stainless steel probe, 12.7mm x 75mm
Core / Strands	7/0.2mm
Cable Termination	Bare Tails
Gland	M20 cable entry gland

Max. Temperature	+150Å°C
Min. Temperature	-50Å°C
IP Rating	67
Accuracy	+/- 0.2K (over 0-70Å°C)
Number of Wires	2