

Prise miniature de connecteur de thermocouple fil rapide AM-E-FQ Type E ANSI



Miniature in-line Socket - Quick Wire
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniature Quick Wire Connecteur Thermocouple Socket AM-E-FQ Type E ANSI

Connexion jab-in facile, il suffit d'enfoncer le fil et de serrer la vis. Les contacts sont polarisés pour éviter une polarisation incorrecte.

Specifications**Specifications**

Product Code	XE-9042-001
General Description	Easy jab-in connection, just push in wire & tighten screw.
Thermocouple Type	E ANSI
Connector Type	Quick Wire Socket
Connector Size	Miniature
Colour	Violet
Max. Temperature	220°C
+Positive Contact	Nickel Chromium
-Negative Contact	Copper Nickel
Standards Met	BS EN 50212:1997. RoHS