



## Connecteur de thermocouple miniature rectangulaire AM-K-FF Type K ANSI



**Miniature Fascia Socket (with nylon mounting clip)**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**connecteur de thermocouple rectangulaire miniature prise de bord AM-K-FF Type K ANSI**

Clip en nylon pour une fixation rapide et simple au panneau. Toutes les douilles de montage du fascia sont polarises au moyen de broches de taille ingale

**Spcifications****Specifications**

<b>Product Code</b>	XE-1137-001
<b>General Description</b>	All fascia mounting sockets are polarised by means of unequally sized pins
<b>Thermocouple Type</b>	K ANSI
<b>Connector Type</b>	Fascia mounting socket
<b>Connector Size</b>	Miniature
<b>Colour</b>	Yellow
<b>Dimensions</b>	See Schematic Drawing (Product Literature Tab)
<b>Max. Temperature</b>	120°C
<b>+Positive Contact</b>	Nickel Chromium
<b>-Negative Contact</b>	Nickel Alloy
<b>Standards Met</b>	BS EN 50212:1997. RoHS