


Miniature Rectangulaire Fascia Socket Thermocouple Connector AM-N-FF Type N ANSI



Miniature Fascia Socket (with nylon mounting clip)
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniature Rectangulaire Thermocouple Connecteur Fascia Socket AM-N-FF Type N ANSI

Clip en nylon pour une fixation rapide et simple au panneau. Toutes les prises de montage fascia sont polarisées au moyen de broches de taille inégale

Spécifications**Specifications**

Product Code	XE-1277-001
General Description	All fascia mounting sockets are polarised by means of unequally sized pins
Thermocouple Type	N ANSI
Connector Type	Fascia mounting socket
Connector Size	Miniature
Colour	Orange
Max. Temperature	120°C
+Positive Contact	Nickel / Chromium / Silicon
-Negative Contact	Nickel Silicon
Standards Met	BS EN 50212:1997. RoHS