


# Miniature Panel Mount Thermocouple Connector with Stainless Steel Bracket AM-R/S-SSPF Type R/S ANSI



**Miniature in-line Socket with stainless steel mounting bracket**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Miniature Panel Mount Thermocouple Connector with Stainless Steel Bracket AM-R/S-SSPF Type R/S ANSI**

Stainless steel bracket simply fitted to connector. Fully enclosed terminals with cable port seal

**Specifications****Specifications**

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | XE-1318-001                                   |
| <b>General Description</b> | Fully enclosed terminals with cable port seal |
| <b>Thermocouple Type</b>   | R/S ANSI                                      |
| <b>Connector Type</b>      | Socket with SS Bracket                        |
| <b>Connector Size</b>      | Miniature                                     |
| <b>Colour</b>              | Green   |
| <b>Max. Temperature</b>    | 220°C   |
| <b>+Positive Contact</b>   | Copper  |
| <b>-Negative Contact</b>   | Copper Alloy                                  |
| <b>Standards Met</b>       | BS EN 50212:1997. RoHS                        |