

## Connecteur de thermocouple standard AS-K-M Type K ANSI



**Standard in-line Plug**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Connecteur de thermocouple standard AS-K-M Type K ANSI**

Les contacts sont fabriqués à partir de véritables alliages de thermocouple. Convient aux applications industrielles

**Spifications****Specifications**

<b>Product Code</b>	XE-1430-001
<b>General Description</b>	Contacts are polarised to prevent incorrect connection
<b>Thermocouple Type</b>	K ANSI
<b>Connector Type</b>	Plug
<b>Connector Size</b>	Standard
<b>Colour</b>	Yellow
<b>Dimensions</b>	See Schematic Drawing (Product Literature Tab)
<b>Max. Temperature</b>	220°C
<b>Pins</b>	Solid Round Pins
<b>+Positive Contact</b>	Nickel Chromium
<b>-Negative Contact</b>	Nickel Alloy
<b>Standards Met</b>	BS EN 50212:1997. RoHS