


Prise de connecteur de thermocouple de fascia rectangulaire standard AS-K-FF Type K ANSI



Standard Fascia Socket (with nylon mounting clip)
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Connecteur de thermocouple rectangulaire standard Femelle de faade AS-K-FF Type K ANSI

Clip en nylon pour une fixation simple et rapide au panneau. Toutes les douilles de montage de faade sont polarises au moyen de broches de taille ingale

Spcifications**Specifications**

Product Code	XE-1437-001
General Description	Accepts standard round pin plug connector. Nylon clip for quick, simple fixing to panel
Thermocouple Type	K ANSI
Connector Type	Fascia Mounting Socket
Connector Size	Standard
Colour	Yellow
Max. Temperature	120°C
+Positive Contact	Nickel Chromium
-Negative Contact	Nickel Alloy
Standards Met	BS EN 50212:1997. RoHS
Series	Standard Connector