


Socket standard rectangulaire de connecteur de fascia thermocouple AS-N-FF Type N ANSI



Standard Fascia Socket (with nylon mounting clip)
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Connecteur thermocouple rectangulaire standard Fascia Socket AS-N-FF Type N ANSI

Clip en nylon pour une fixation rapide et simple au panneau. Toutes les prises de montage fascia sont polarisées au moyen de broches de taille inégale

Spécifications**Specifications**

Product Code	XE-1577-001
General Description	Accepts standard round pin plug connector. Nylon clip for quick, simple fixing to panel
Thermocouple Type	N ANSI
Connector Type	Fascia Mounting Socket
Connector Size	Standard
Colour	Orange
Max. Temperature	120°C
+Positive Contact	Nickel / Chromium / Silicon
-Negative Contact	Nickel Silicon
Standards Met	BS EN 50212:1997. RoHS