


## Socket standard rectangulaire de connecteur de fascia thermocouple AS-R/S-FF Type R/S ANSI



**Standard Fascia Socket (with nylon mounting clip)**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Connecteur thermocouple rectangulaire standard Fascia Socket AS-R/S-FF Type R/S ANSI**

Clip en nylon pour une fixation rapide et simple au panneau. Toutes les prises de montage fascia sont polarisées au moyen de broches de taille inégale

**Spécifications****Specifications**

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | XE-1617-001   |
| <b>General Description</b> | Accepts standard round pin plug connector. Nylon clip for quick, simple fixing to panel |
| <b>Thermocouple Type</b>   | R/S ANSI  |
| <b>Connector Type</b>      | Fascia Mounting Socket  |
| <b>Connector Size</b>      | Standard  |
| <b>Colour</b>              | Green   |
| <b>Max. Temperature</b>    | 120°C   |
| <b>+Positive Contact</b>   | Copper  |
| <b>-Negative Contact</b>   | Copper Alloy  |
| <b>Standards Met</b>       | BS EN 50212:1997. RoHS  |