



# Miniature Quick Wire Thermocouple Connector Plug JM-B-MQ Type B JIS



**Miniature in-line Plug - Quick Wire**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Miniature Thermocouple Connecteur Quick Wire Plug JM-B-MQ Type B JIS**

Connecteur miniature avec broches a bouchon plat. Connexion jab-in facile pour la résiliation rapide. Il suffit de pousser le fil et serrer la vis.

**Spécifications****Specifications**

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | XE-9222-001                                   |
| <b>General Description</b> | Easy jab-in connection for rapid termination. |
| <b>Thermocouple Type</b>   | B JIS   |
| <b>Connector Type</b>      | Quick Wire Plug                               |
| <b>Connector Size</b>      | Miniature                                     |
| <b>Colour</b>              | Grey  |
| <b>Max. Temperature</b>    | 220°C   |
| <b>Pins</b>                | Flat Pins                                     |
| <b>+Positive Contact</b>   | Copper  |
| <b>-Negative Contact</b>   | Copper  |
| <b>Standards Met</b>       | BS EN 50212:1997. RoHS                        |