

## Prise miniature de connecteur de thermocouple fil rapide JM-T-FQ Type T JIS



**Miniature in-line Socket - Quick Wire**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Miniature Quick Wire Connecteur Thermocouple Socket JM-T-FQ Type T JIS**

Connexion jab-in facile, il suffit d'enfoncer le fil et de serrer la vis. Les contacts sont polarisés pour éviter une polarisation incorrecte.

**Specifications****Specifications**

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | XE-9263-001  |
| <b>General Description</b> | Easy jab-in connection, just push in wire & tighten screw. |
| <b>Thermocouple Type</b>   | T JIS  |
| <b>Connector Type</b>      | Quick Wire Socket  |
| <b>Connector Size</b>      | Miniature  |
| <b>Colour</b>              | Brown  |
| <b>Max. Temperature</b>    | 220°C  |
| <b>+Positive Contact</b>   | Copper   |
| <b>-Negative Contact</b>   | Copper Nickel  |
| <b>Standards Met</b>       | BS EN 50212:1997. RoHS                                     |