

# Miniature Rectangulaire Fascia Socket Thermocouple Connecteur JM-R/S-FF Type R/S JIS



**Miniature Fascia Socket (with nylon mounting clip)**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Miniature Rectangulaire Thermocouple Connecteur Fascia Socket JM-R/S-FF Type R/S JIS**

Clip en nylon pour une fixation rapide et simple au panneau. Toutes les prises de montage fascia sont polarisées au moyen de broches de taille inégale

**Spécifications****Specifications**

<b>Product Code</b>	XE-9197-001
<b>General Description</b>	All fascia mounting sockets are polarised by means of unequally sized pins
<b>Thermocouple Type</b>	R/S JIS
<b>Connector Type</b>	Fascia mounting socket
<b>Connector Size</b>	Miniature
<b>Colour</b>	Black
<b>Max. Temperature</b>	120°C
<b>+Positive Contact</b>	Copper
<b>-Negative Contact</b>	Copper Alloy
<b>Standards Met</b>	BS EN 50212:1997. RoHS