

Miniature Panel Mount Thermocouple Connector avec Stainless Steel Bracket JM-R/S-SSPF Type R/S JIS



Miniature in-line Socket with stainless steel mounting bracket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniature Panel Mount Thermocouple Connector avec Stainless Steel Bracket JM-R/S-SSPF Type R/S JIS

Support en acier inoxydable simplement monté sur connecteur. Terminaux entièrement fermés avec joint de câble

Spécifications**Specifications**

Product Code	XE-9198-001
General Description	Fully enclosed terminals with cable port seal
Thermocouple Type	R/S JIS
Connector Type	Socket with SS Bracket
Connector Size	Miniature
Colour	Black
Max. Temperature	220°C
+Positive Contact	Copper
-Negative Contact	Copper Alloy
Standards Met	BS EN 50212:1997. RoHS