



# Miniature Panel Mount Thermocouple Connector avec Stainless Steel Bracket JM-B-SSPF Type B JIS



**Miniature in-line Socket with stainless steel mounting bracket**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Miniature Panel Mount Thermocouple Connector avec Stainless Steel Bracket JM-B-SSPF Type B JIS**

Support en acier inoxydable simplement monté sur connecteur. Terminaux entièrement fermés avec joint de câble

**Spécifications****Specifications**

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | XE-9218-001                                   |
| <b>General Description</b> | Fully enclosed terminals with cable port seal |
| <b>Thermocouple Type</b>   | B JIS   |
| <b>Connector Type</b>      | Socket with SS Bracket                        |
| <b>Connector Size</b>      | Miniature                                     |
| <b>Colour</b>              | Grey  |
| <b>Max. Temperature</b>    | 220°C   |
| <b>+Positive Contact</b>   | Copper  |
| <b>-Negative Contact</b>   | Copper  |
| <b>Standards Met</b>       | BS EN 50212:1997. RoHS                        |