



Standard Thermocouple Connector Plug JS-R/S-M Type R/S JIS



Standard in-line Plug
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Standard Thermocouple Connector Plug JS-R/S-M Type R/S JIS

Les contacts sont faits a partir de véritables alliages Thermocouple. Convient aux applications industrielles

Spécifications**Specifications**

| | |
|----------------------------|--|
| Product Code | XE-9310-001 |
| General Description | Contacts are polarised to prevent incorrect connection |
| Thermocouple Type | R/S JIS |
| Connector Type | Plug |
| Connector Size | Standard |
| Colour | Black |
| Max. Temperature | 220°C |
| Pins | Solid Round Pins |
| +Positive Contact | Copper |
| -Negative Contact | Copper Alloy |
| Standards Met | BS EN 50212:1997. RoHS |