


# Standard Quick Wire Panel Mount Thermocouple Connector avec Bracket en acier inoxydable JS-J-SSPFQ Type J JIS



**Standard Socket - Quick Wire with stainless steel mounting bracket**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Standard Quick Wire Panel Mount Thermocouple Connector avec Bracket en acier inoxydable JS-J-SSPFQ  
Type J JIS**

Destiné à être utilisé dans des procédés plus industriels, mais peut être utilisé dans toutes les applications comme souhaité. Connexion facile 'jab-in' pour la terminaison rapide, il suffit de pousser le fil et serrer la vis

**Spécifications****Specifications**

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | XE-9299-001   |
| <b>General Description</b> | Easy 'jab-in' connection for rapid termination, just push wire in and tighten screw |
| <b>Thermocouple Type</b>   | J JIS   |
| <b>Connector Type</b>      | Quick Wire Socket with SS Bracket   |
| <b>Connector Size</b>      | Standard  |
| <b>Colour</b>              | Yellow  |
| <b>Dimensions</b>          | See Schematic Drawing (Product Literature Tab)                                      |
| <b>Max. Temperature</b>    | 220°C   |
| <b>+Positive Contact</b>   | Iron  |
| <b>-Negative Contact</b>   | Copper Nickel   |
| <b>Standards Met</b>       | BS EN 50212:1997. RoHS  |