

Barrier Terminal Block Thermocouple Connector JS-R/S-BTS Type R/S JIS



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Barrier Terminal Block Thermocouple Connecteur JS-R/S-BTS Type R/S JIS

Fourni en paires simples et imbriquées. Connexions dans le matériel thermocouple

Voir l'onglet spécifications pour tous les détails .

Specifications

Product Code	XE-1963-001
General Description	Supplied in single, interlocking pairs. Connections in thermocouple material
Thermocouple Type	R/S JIS
Connector Type	Barrier Terminal Blocks
Colour	Black
Dimensions	See Schematic Drawing (Product Literature Tab)
Max. Temperature	220°C
+Positive Contact	Copper
-Negative Contact	Copper Alloy
Standards Met	BS EN 50212:1997. RoHS