

## Barrier Terminal Block Thermocouple Connector JS-R/S-BTS Type R/S JIS



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Barrier Terminal Block Thermocouple Connecteur JS-R/S-BTS Type R/S JIS**

Fourni en paires simples et imbriquées. Connexions dans le matériel thermocouple

**Voir l'onglet spécifications pour tous les détails .**

## Specifications

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | XE-1963-001  |
| <b>General Description</b> | Supplied in single, interlocking pairs. Connections in thermocouple material |
| <b>Thermocouple Type</b>   | R/S JIS  |
| <b>Connector Type</b>      | Barrier Terminal Blocks  |
| <b>Colour</b>              | Black  |
| <b>Dimensions</b>          | See Schematic Drawing (Product Literature Tab)                               |
| <b>Max. Temperature</b>    | 220°C  |
| <b>+Positive Contact</b>   | Copper   |
| <b>-Negative Contact</b>   | Copper Alloy   |
| <b>Standards Met</b>       | BS EN 50212:1997. RoHS   |