


illet de dcharge de traction



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Passe-cble de dcharge de traction pour connecteurs miniatures

Temprature nominale : -30C +140C

Matriau : lastomre thermoplastique (TPE)

Specifications

| | |
|---------------------|-------------|
| Product Code | XE-1882-001 |
|---------------------|-------------|